Email completed form to Quotes@RehabMart.com Phone: 1-800-827-8283

## Convaid Metro - HCPCS E1037

				Date:	
Contact:	Mark For:				
Bill To:	Ship To:				
Address:	Address:				
City/State/Zip:					
Phone:	City/State/Zip:				
Email:	Phone:				
Metro Standard Features:	•				
To 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (ME12/14/16) Quick release Rolko rear wheels (ME12/14/16)  R" Front, 11" Rear Solid Tire (ME18)	<ul> <li>30° Fixed Tilt</li> <li>Silver Powder Coated Frame</li> <li>2-Pt Positioning Belt</li> <li>Made in the USA</li> </ul>				
Upholstery Colors  ☐ Cherry Red		ME12 C903850	ME14 C903850	ME16 C903850	ME18 C903850
☐ Royal Blue		C903855	C903855	C903855	C903855
Tire Options  ☐ 12.5" Rear Solid Knobby Tire  ☐ 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire  ☐ 8"x2" Front Ribbed, 12.5" Rear Solid Knobby Tire  ☐ 8"x2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire	<u>HCPCS</u>	<b>ME12</b> C904117 C904118	ME14 C904117 C904118	ME16 C904117 C904118	ME18  C903084  C903083
Transit Required Items	<u>HCPCS</u>	ME12	ME14	ME16	ME18
<ul> <li>□ Transit Option (Manufacturer Installation Required)</li> <li>□ Headrest Extension</li> <li>□ 3-Pt Positioning Belt</li> <li>Foot Positioners - Select Size Below:</li> <li>□ Foot Positioners - Small 7" - 10" (pair)</li> <li>□ Foot Positioners - Medium 8" - 11" (pair)</li> <li>□ Foot Positioners - Large 9" - 12" (pair)</li> </ul>	E0966 E0978	C901403 C903263 C904017 C903415 C903416 C903417	C901401 C903264 C904017 C903415 C903416 C903417	C901336 C903265 C904017 C903415 C903416 C903417	C903128 C903266 C904017 C903415 C903416 C903417
Transit Ontions	HCPCS	ME12	ME14	ME16	ME18
Transit Options  ☐ Q'Straint Transit Lap Belt ☐ Sure-Lok Transit Lap Belt	E0978 E0978	C903094 C903232	C903094 C903232	C903094 C903232	C903094 C903232
Accessories	HCPCS	ME12	ME14	ME16	ME18
☐ Wheelchair Travel Bag ☐ Under Seat Storage Basket - Fixed Tilt		C902022 C904046	C902022 C904047	C902022 C904047	C902022 C904048

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The Metro Transit option includes four manufacturer installed wheelchair transportation anchors. The Metro transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are ME12: 66 lbs./30kg; ME14: 100 lbs./45.5kg/ME16 and ME18: 170 lbs./77 kg. The Convaid Metro was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is being, or can be, made.



1-800-827-8283 3651 Mars Hill Rd, Bldg 2400 Watkinsville, GA 30677-6042 USA

Page 1 of 1 12/20/2022