

# Odyssey Chair

The Odyssey Chair is a multifunctioning and highly adjustable chair designed to meet the needs of physically challenged children not only today but in the future as well. There are two sizes available: Early Childhood and Secondary. As an added benefit there is an upgrade package available that allows the Early Childhood Chair to grow to a Secondary size chair without purchasing a whole new chair!

The modular concept provides multiple options and adjustments with tremendous "growth" potential. It is an ideal chair for evaluation, treatment, and seating intervention in the home, school, or clinic setting. The easy adjustability of the chair can meet the changing needs of the user over the years.

The sides of the chair are constructed of laminated phenolic to withstand years of use. The wooden seat is height adjustable within the chair to allow it to be positioned low for self transfers, at a mid-range height to use at a school desk or play table, or up high as a feeding chair at the kitchen or dining room table. The seat is also depth and angle adjustable for individualized therapeutic pelvic and lower extremity positioning. The height and depth adjustable back can also be angled 30° back in space.





BASE UNIT



Addie, age 5, is able to play or eat at the family table in the Odyssey Chair (THA-OC-150).

## FIT RIGHT GUIDE

See center insert page L



### Odyssey Chair, Choose Size

#### EARLY CHILDHOOD

- THA-OC-150 Odyssey Chair
- Seat depth (B1): 7" - 14"
- Seat width (G): 12½" (6" to 12½")\*
- Seat height to footrest (A): 5" to 14"
- Seat height to floor: 11" - 19"\*\*\*
- Back height (D): 15"
- THA-OC-150AB Abduction Block
- THA-OC-150AR Armrests (pair)
- Armrest height (C1): 6" - 10"
- THA-OC-150BK Butterfly Kit
- THA-OC-150FS Foot Sandals (pair)
- Inside dimensions: 3"W x 5"L
- Outside dimensions: 4¼"W x 5½"L
- THA-OC-150HR Headrest
- Back height (E): 21" - 26"
- THA-OC-150HS Heel Stop
- THA-OC-150PS Seat and Back Pad Set
- THA-OC-150PSW Seat (full width) and Back Pad Set
- THA-OC-150SB Support Blocks\* 5"H x 4"W (pair)
- THA-OC-150SK Support Kit
- (includes FS, PS, and 2 pair SB)
- THA-OC-150SW Head Support Wedges 6"H x 6"W (pair)
- THA-OC-150T Wood Tray\*\*\*
- Tray dimensions: 20"W x 18"D
- Tray cutout: 14"W x 4"D
- THA-OC-150TC Rounded Clear Polycarbonate Tray\*\*\*
- Tray dimensions: 26"W x 21"D
- Tray cutout: 13½"W x 6"D
- THA-OC-150WK Locking Wheel Kit
- Height: 3"H
- THA-OC-150WK-FL Full-Lock Wheel Kit
- Height: 4"H

\*Note: Seat width reduction is achieved by use of Support Blocks.

\*\*Note: Height measured without Wheel Kit.

\*\*\*Note: Armrests must be ordered to attach the Tray to the chair.

THA-OC-UP Upgrade Kit   
(Expand THA-OC-150 to THA-OC-250)

**SECONDARY**

- THA-OC-250 Odyssey Chair □
- Seat depth (B1): 11" - 16" □
- Seat width (G): 15½" (8" to 15½")\* □
- Seat height to footrest (A): 5" to 18" □
- Seat height to floor: 11" - 19"\*\*\* □
- Back height (D): 19" □
- THA-OC-250AB Abduction Block □
- THA-OC-250AR Armrests (pair) □
- Armrest height (C1): 6" - 10" □
- THA-OC-250BK Butterfly Kit □
- THA-OC-250FS Foot Sandals (pair) □
- Inside dimensions: 4"W x 8"L □
- Outside dimensions: 5¼"W x 8½"L □
- THA-OC-250HR Headrest □
- Back height (E): 25" - 30" □
- THA-OC-250HS Heel Stop □
- THA-OC-250PS Seat and Back Pad Set □
- THA-OC-250PSW Seat (full width) and Back Pad Set □
- THA-OC-250SB Support Blocks\* 5"H x 6"W (pair) □
- THA-OC-250SK Support Kit □
- (includes FS, PS, and 2 pair SB) □
- THA-OC-250SW Head Support Wedges 6"H x 6"W (pair) □
- THA-OC-250T Wood Tray\*\*\* □
- Tray dimensions: 25"W x 21"D □
- Tray cutout: 17"W x 4"D □
- THA-OC-250TC Rounded Clear Polycarbonate Tray\*\*\* □
- Tray dimensions: 26"W x 21"D □
- Tray cutout: 13½"W x 6"D □
- THA-OC-250WK Locking Wheel Kit □
- Height: 3"H □
- THA-OC-250WK-FL Full-Lock Wheel Kit □
- Height: 4"H □

\*Note: Seat width reduction is achieved by use of Support Blocks.  
 \*\*Note: Height measured without Wheel Kit.  
 \*\*\*Note: Armrests must be ordered to attach the Tray to the chair.



*The stability of the Odyssey Chair and the strength of its footplate allow Addie to climb into the chair independently.*



**OPTIONS**

- A.** The cost effective **Upgrade Kit** converts the Early Childhood Chair to a Secondary Chair using the components you already have. It comes complete with a new wood seat, back, and footboard, plus the hardware necessary to increase the width of the chair from 12½" – 15½" wide.
- B.** A **Heel Stop**, complete dual ankle straps and curved side release buckles, is available to help keep the feet secure on the footplate. It is ideal for users who will be standing on the foot board and transferring in and out of the chair.
- C.** Two **Wheel Kits**, both complete with mounting plates and four independent locking wheels, increase the portability of the chair from room to room, allowing the child to remain seated and well positioned for longer periods of time. The **Locking Wheel Kit's** toggle lock prevents the wheels from rolling movement but not swivel. The **Full-Lock Wheel Kit's** wheels lock, stopping both rolling and swivel movement, making it a more stable choice for independent or assisted transfers.
- D.** See pages 20-21 for these options shown: Abduction Block, Armrests, Butterfly Kit, Foot Sandals, Headrest, Seat and Back Pad Set, Support Blocks (4), Head Support Wedges.