

**Defender Reha US Car Seat** – Jan. 2024 US Retail prices

Date	Buyer		PO#		Mark for
Address					
Contact		Phone		Fax	

## Step 1 - Seat



HCPCS

Description

Seat and back with height-adjustable headrest, LATCH equipped, 5-point safety/ 99 positioning belt with shoulder pads, chest clip, small crotch pad, set of removable head, trunk and hip supports, cup holder, color black-charcoal, 2 years warranty











Step 2 - General Options

•		
General Accessories	HCPCS	Unit
Footrest, short (depth 4.7")* (angle adjust.)	E1399	Each
Footrest, long (depth 7.1")* (angle adjust.)	E1399	Each
Footrest adapter + Base plate 89036*	E1399	Each
Seat wedge, below, black	E1399	Each
Crotch pad, large, black**	E1399	Each
Anterior support, table, black	E1399	Each
Abduction block, padded, black	E1399	Each
Footrest pad (for 806), height 4", black	E1399	Each
Front stand	E1399	Each

Technical Specs Defender Reha

Child weight	see user manual		
Seat weight	22 lbs		
Seat depth	13.5"		
Seat width	8.5-12" ***		
Seat to head	23.5-30.5" ****		
Shoulder height	10-19.5" *****		
Lower leg length	(7")# 11-17"		

#7"minimum with footrest pad 834

This is a Universal Restraint System for disabled children. It is approved in accordance with Federal Motor Vehicle Safety Standard 213 for general use in vehicles and it will fit most but not all cars. A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a Universal Restraint System for this age group. Only suitable for use in **vehicles with 3-point** (lap & shoulder) seat belts. This is a seat for children weighing between 22 and 65 lbs., whose height is between 27 and 57 inches and who are over 12 month old in **Harness Mode** or 40 – 110 lbs., and whose height is between 34 and 57 inches in **Booster Mode**. Suitable only for **forward-facing** installation.

- \* The footrest adapter 821 & Base plate 89036  $\underline{\text{must}}$  be ordered to use footrest
- \*\* We recommend using larger crotch pad only in booster mode
- \*\*\* To reach max. seat width remove additional hip guide pads from additional seat pad
- \*\*\*\* To reach max. back height remove additional seat pad and foam pad under seat cushion
- \*\*\*\*\* To reach max. shoulder height remove additional seat pad and foam pad under seat cushion







Table



Seat wedge below



Abduction block



Crotch pad