

Plastic Cylinder Wrenches (E-Keys)

EK-1 25/pk
Also known as "E Keys".
Used on Post Valve
Oxygen Cylinder
to Open and Close.

COMPARE TO:

	Hudson®	
EK-1	50893	




E-Key

Oxygen Regulator Connector Swivel

XM-1 (XMAS Trees), 25/pk
For Use on
Concentrator or
Regulator with DISS
Fitting in Order to
Attach Cannula.

COMPARE TO:

	Hudson® Salter®	
XM-1	2555	1205



XM-1

Cylinder Wrench

18501 For Use with Small Cylinders
used in Most Portable Oxygen
Set Ups, 10/pkg

FEATURES:

- Has place to store Viton rubber/brass sealing washer.
- Used to open post valves.
- Can be used to assist in install or removal of both large wing nut and small t-bar regulators.



18501

Concentrator Water Trap

CON 900 50/bx

FEATURES:

- Easy inline downstream installation with small bore tubing. Can be cleaned with detergent or in a vinegar bath.

COMPARE TO:

	Hudson® Salter®	
CON 900	1679	7000



CON 900

Tubing Extension Connector


CON 400 Swivel, 50/bg
CON 500 Straight, 50/bg
CON 550 Ridged, 50/bg
CON 700 Y Connector, 10/bg

FEATURES:

- Rigid plastic tubing connector allows extension tubing to be connected to present cannula for greater length(s).

COMPARE TO:

	Hudson® Salter®	
CON 400	N/A	1220
CON 500	N/A	2005
CON 550	1420	1215
XM-1	2555	1205



CON 400

CON 500

CON 550

CON 700