




PHYSIO MAGNETO SUPER TRANSDUCTION
PLUS LASER THERAPY

FIELD INTENSITY UP 30%
DIODO LASER POWER MAX 5600MW



PMST NEO+

PMST + LASER



100% Safe and Painless and Effective :

Security and stability of home laser as well as the exact effect of medical laser. The treatment process is painless and harmless, non-invasive, no trauma, 100% safe. Laser therapy delivers non-thermal (non-heat) photons of light to the body to help heal injured cells.

Cold DIODO Laser Therapy :

A perfect combination of medical laser **808nm** and home laser **650nm** will have better curative effect for body pain relief. this Laser which is the best for the treatment of humans and animals.

Cold DIODO Laser Therapy Advantages :

Proven results for pain relief, inflammation, soft tissue, and skin injury therapy. Natural healing, non-invasive, no side effects and no cross infection.

Multifunctional, Portable,
Max Power Of **5600MW**.

PHYSIO MAGNETO SUPER TRANSDUCTION

PMST NEO+

FIELD INTENSITY UP







The PMST NEO + opens up new possibilities in regeneration and rehabilitation. The areas of application include diseases of the musculoskeletal system, joints pain, sports injuries, muscle strain, fracture healing etc.

PMST differs from general forms of magnetic field therapy or PEMF due to the high oscillation frequency, This property enables a high penetration depth and a large range of indications.

- PAIN FREE PHYSIOTHERAPY
- HIGH INTENSITY MAGNETIC FIELD STRENGTH
- COMFORTABLE AND NO NEED TO UNDRESS
- HANDS FREE WORK WITH HELP OF HOLDING ARM
- COMBINE WITH DIODO LASER THERAPY
- HIGH EFFECTIVE OF WATER COOLING SYSTEM

PMST + LASER



INTEGRATED APPLICATOR

FOR BETTER PHYSIOTHERAPY EFFECTS



The PMST NEO+ features unique applicator design. The ring type electromagnetic coil applicator connect with LASER applicator by the special design connector. It is the only one of its kind in world physiotherapy field , can transduce magnetic pulse deep into body tissue ,at the same time, DIODO LASER focused on the same treatment area. The two technologies perfectly combine together for better therapeutic effects.



PMST + LASER

Pulsed Electromagnetic Field Therapy

PMST different with PEMF, it is a ring type coil, cover bigger area and fit to joints part. High speed oscillation for deeper penetration.



PMST NEO+ LASER

PHYSIO MAGNETO SUPER TRANSDUCTION PLUS LASER THERAPY OF 5600MW

- PHYSIO MAGNETO
- Frequency up to 3000Hz
- Intensity up to 5 Tesla
- Unique integrated applicator
- Flexible holding arm




Intuitive operation with
10.4 inch touch screen



3 joints holding arm
Easy applicator positioning

PMST + LASER

UNIQUE RING TYPE FIELD ▶



◀ DIODO COLD LASER THERAPY

TECHNICAL SPECIFICATIONS OF PMST NEO+

Field strength at the coil	5T
Voltage	100-240V
Output power	1500W
Oscillation frequency in ST mode	1000-3000Hz
Oscillation frequency in MT mode	1-100Hz
Water cooling system	Water 2.5L
Diodo cold laser wavelength	808NM / 650NM
Package	Alu box inside and carton box outside
Package size	66*60*50KG
Gross weight	44K