

Child transfer to Churchill Base

- Easy transfer of child to vehicle
- Super low profile base
- Swing away abductor
- Shoulder, hip and thigh width unlimited
- 15" or 18" seat depth available



Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Churchill Belt Guides



- Optional belt guides to position vehicle seat belt low on upper thighs
- For larger child or adult belt guides can be moved to side if not needed – this will allow more hip room

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Churchill Vest /Shoulder Belt Guide



Easy to use shoulder belt guides to keep shoulder belt off the neck and positioned properly over the chest



Easy to install and adjust vest for upper body control
(Two vest sizes available)

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Hip Inflection Wedge



Hip inflection wedge



Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Churchill Pommels



- Standard pommel for keeping legs apart and in position



- Comfort Pommel to keep legs apart and position legs into a more comfortable position

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Leg Straps to Control Slouching



- Adjustable leg straps (standard)
- Leg straps keep child / adult hips in position – prevent slouching
- Can be used with or without the pommel

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

A-ductor Strap

- Optional A-ductor strap to position and control leg movement
- A-ductor strap can be used with or without the pommel



Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Churchill Head Rest



- Easy to install headrest mounts to vehicle head rest
- Provides lateral head support and sleeping support aid
- When vehicle does not have a headrest, the Churchill headrest attaches to seat mounting strap

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Churchill EZ-Up Caps



Available in three sizes

- EZ-Up Cap adjustable to fit the child's head
- EZ-Up Cap attaches to velcro located on headrest to keep head from falling forward

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™

Latch Bracket Attachment Kit



The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH.

The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket Kit will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the Attachment Bracket Kit can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket Kit cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.

Merritt Manufacturing Inc.

Listening to Therapists; Working for Families™