Protekt® Positioning Devices

ROLL CONTROL BOLSTER

- Alternative to side rails.
- Prevents the patient from falling out of bed.
- Easily attaches to bed with quick-release buckles and can be positioned to accommodate individual patient needs.
- Provides secure positioning for side-lying patients.
- Not considered a restraint.

| 58500 | Protekt® Roll Control Bolster | Pair 34"x7"x8" |
|-------|-------------------------------|----------------|
| 36300 | FIGURES ROIL CONTROL BOISTEI | raii 34 X/ X0 |



pain.

- Improves positioning & circulation.
- High density polyurethane foam construction.
- Cover Options: Removable fluid resistant nylon or vinyl.

| 50010 | Protekt® Leg Elevator Wedge | 20"x24"x6" |
|-------|-----------------------------|-------------|
| 50011 | Protekt® Leg Elevator Wedge | 20"x24"x8" |
| 50012 | Protekt® Leg Elevator Wedge | 20"x24"x10" |

BED WEDGE

- Designed to be used as a multipurpose wedge and can be placed in a variety of support positions.
 Ideal for head, leg or foot elevation.
- The gradual slope helps ease respiratory problems while reducing neck, back and shoulder pain.
- Improves positioning and circulation.
- High density polyurethane foam construction.
- Cover: Breathable polyester/cotton fabric.

| 50004 | Protekt® Foam Bed Wedge 7" | 23"x23"x7" |
|-------|-----------------------------|-------------|
| 50005 | Protekt® Foam Bed Wedge 10" | 23"x23"x10" |
| 50006 | Protekt® Foam Bed Wedge 12" | 23"x23"x12" |



