

PORTABLE OZONE CPAP/BiPAP SANITIZER

Purify Q

ON CHARGE

DISINFECT WITH JUST THE TOUCH OF A BUTTON

#190-6000

The simple & effective way to clean your CPAP/BiPAP mask and equipment

Powerful, safe hospital grade sanitation cycle utilizes ozone to kill 99.99% of germs and bacteria in just 35 minutes! Lightweight & compact - ideal for use at home and away.

No extra supplies required - sanitizes and stores in the included cleaning/travel case!

Don't take a chance with your health! Sanitize with \mathcal{P} urify \boldsymbol{Q}







SPECIFICATIONS:

Size/Dimensions: 4.75" x 3.75" x 2.12" Weight: 0.5 lb Cleaning Technology: Ozone (O3) Power Source: 4W, Rechargeable Lithium-ion battery Recommended Frequency of Use:

- Three times per week: mask sanitation
- Once per week: humidification chamber, tubing
- As needed: nebulizer mask, tubing, breast pump accessories, other medical devices

Includes: Purify O3 unit, dual purpose sanitizing and travel case, tubing adapter, USB cable and charging adapter, bilingual instruction manual (English and Spanish) Warranty: 2 years - register at *respondo2.com/warranty* FDA Registered

SIMPLE TO USE:

- Empty humidifier chamber and wipe mask/tubing to remove visible residue, attach and or place Purify O3 in sanitizing case or closed container.
- Closed container disinfection process ensures safe and effective equipment sanitation at a higher efficacy than other cleaning methods.
- Push to start Easy one-touch control. Automatically shuts off after 35 minute cleaning cycle is complete. (Note: Allow ozone to dissipate for additional 90 minutes before use.)
- Recharges quickly via wall outlet or USB port.
- No maintenance, filters, or extra cleaning solutions required.

PURIFY O3 CPAP/BiPAP SANITIZER KIT

#190-6000 Purify O3 Sanitizer Kit

REPLACEMENT COMPONENT PARTS

#190-6002Dual Purpose Sanitizing/Travel Case#190-6003USB Cable#190-6004Charging Adapter#190-6005Tubing Adapter

