

Midline

Positioning Wheelchair

POSITIONING
BRODA[®]




Full tilt and recline functionality provides superior positioning



Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt-in-space enhances comfort



Comfort Tension Seating[®] increases support as well as sitting tolerance

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Allows for easy maneuverability and safe patient adjustments
- Reduces risk of injury and falls with supportive stable positioning

*Enhanced positioning
with refined features*

SPECIFICATIONS

<i>Degrees of Tilt</i>	<i>Up to 38°</i>
<i>Recline</i>	<i>Up to 90° (Trendelenburg capable)</i>
<i>Weight Capacity</i>	<i>up to 350 lbs. (114 kg)</i>
<i>Seat Widths</i>	<i>16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm)</i>
<i>Seat Depths</i>	<i>15.5", 17", 18.5", 20" (39cm, 43cm, 45cm, 58cm)</i>
<i>Floor-to-Seat Heights</i>	<i>Adjustable 17", 18", 19", 20", 21", 22" (43cm - 55cm)</i>
<i>Wheels</i>	<i>5" Front Casters: Total Lock (12.5cm) 5" Rear Casters: Directional or Total Lock (20cm)</i>
<i>Lifetime Warranty Frame</i>	<i>10 years on 16 Gauge tubular steel frame</i>
<i>Pressure-Reducing Seat & Backrest</i>	<i>Integrated with Comfort Tension Seating[®]</i>
<i>Standard Padding</i>	<i>Neckrest, shoulder bolsters, upper and lower side pads, seat and back pad, calf and sole pad</i>