

CLASSIC SERIES  
**McKENZIE ELECTRIC TABLES**

---

- McKenzie Basic
- McKenzie Lift Back
- McKenzie Deluxe

# MCKENZIE BASIC ELECTRIC TABLE

## McKenzie Basic Electric Table

MK-2807

**Lifting Capacity: 400 lbs.**

**Dimensions: 73"L x 28"W x 19-38"H**

With our greatest height range yet, our McKenzie Basic Electric Massage Table raises smoothly and effortlessly from 19 to 38 inches, giving you the most options during treatment. Plus, our rock solid construction—complete with retractable wheels—also ensures that your clients feel incredibly comfortable at any height.

- Height range is 19-38"
- Table weight: 140 lbs.
- Dual-action face rest and pillow
- Standard single foot control pedal
- Whisper quiet lift system
- Retractable wheels
- 3" Deluxe wrap foam system
- Durable precision cut steel base
- Dual end accessory ports (9" on center)
- Table is ADA compliant (consult your tax professional about eligibility for a Section 44 ADA Tax Credit)

- .....
- MK-3007** - 30" wide cushion
  - MK-3307** - 33" wide cushion
  - NTU** - Upgrade to Natural Touch Upholstery
  - NT/TRU** - Upgrade to Natural Touch Upholstery and Touch Response Foam
  - PG** - Prenatal
  - BR** - Breast recesses
  - RC** - Rounded corners
  - SEH-H** - Holes for Side Arm Extensions (Head End)
  - SEH-F** - Holes for Side Arm Extensions (Foot End)
  - XL** - Extra long table (78")
  - 8OC** - Custom configuration 8" on center face rest holes
  - CRFH+** - Crescent face hole with pillow and plug
  - UP DPOS-B** - Upgrade to Deluxe Adjustable face rest base
  - UP PPOS-B** - Upgrade to Pivot Posi-tilt face rest base



### Standard Upholstery Color Choices:



### Upholstery Upgrade Natural Touch Choices:



# MCKENZIE LIFT BACK ELECTRIC TABLE

## McKenzie Lift Back Electric Table

TMK-2807

**Lifting Capacity: 400 lbs.**

**Dimensions: 73"L x 28"W x 19-38"H**

The flexibility of a manual lift back enhances the incredible height range in our McKenzie Lift Back Electric Massage Table, which raises smoothly and effortlessly from 19 to 38 inches. Your clients will feel both comfortable and secure due to rock solid construction—including retractable wheels. Add optional swivel side arms or upgrade to an electric lift back for even greater comfort.

- Height range is 19-38"
- Table weight: 140 lbs.
- Dual-action face rest and pillow
- Side arm extension holes (head end - swivel side arm supports not included)
- Standard single foot control pedal
- Whisper quiet lift system
- Retractable wheels
- 3" Deluxe wrap foam system
- Durable precision cut steel base
- Dual end accessory ports (9" on center)
- Table is ADA compliant (consult your tax professional about eligibility for a Section 44 ADA Tax Credit)

- 
- TMK-3007** - 30" wide cushion
  - TMK-3307** - 33" wide cushion
  - NTU** - Upgrade to Natural Touch Upholstery
  - NT/TRU** - Upgrade to Natural Touch Upholstery and Touch Response Foam
  - PG** - Prenatal (Foot end only)
  - BR** - Breast recesses
  - RC** - Rounded corners
  - 8OC** - Custom configuration 8" on center face rest holes
  - CRFH+** - Crescent face hole with pillow and plug
  - UP DPOS-B** - Upgrade to Deluxe Adjustable face rest base
  - UP PPOS-B** - Upgrade to Pivot Posi-tilt face rest base
  - PTB** - Power tilt back



### Standard Upholstery Color Choices:



### Upholstery Upgrade Natural Touch Choices:



# MCKENZIE DELUXE ELECTRIC TABLE

## McKenzie Deluxe Electric Table

YK-2807

**Lifting Capacity: 400 lbs.**

**Dimensions: 73"L x 28"W x 19-38"H**

Our most adjustable table in the McKenzie line, the McKenzie Deluxe Electric Massage Table combines our incredible height range (19 to 38 inches) with the flexibility of a manual back and knee lift. You can move from massages to facials to pedicures with ease. Your clients will always feel secure and comfortable due to our rock-solid construction—at any height—and retractable wheels. (Doesn't include Crescent Face Pillow)

- Height range is 19-38"
- Table weight: 140 lbs.
- Salon Head Support and Neck Pillow
- Dual action face rest base
- Side arm extension holes and swivel side arm supports (head end)
- Standard single foot control pedal
- Whisper quiet lift system
- Retractable wheels
- 3" Deluxe wrap foam system
- Rounded corners
- Durable precision cut steel base
- Dual end accessory ports (9" on center)
- Table is ADA compliant (consult your tax professional about eligibility for a Section 44 ADA Tax Credit)

.....

- YK-3007** - 30" wide cushion
- YK-3307** - 33" wide cushion
- NTU** - Upgrade to Natural Touch Upholstery
- NT/TRU** - Upgrade to Natural Touch Upholstery and Touch Response Foam
- BR** - Breast recesses
- 8OC** - Custom configuration 8" on center face rest holes
- CRFH+** - Crescent face hole with pillow and plug
- UP DPOS-B** - Upgrade to Deluxe Adjustable face rest base
- UP PPOS-B** - Upgrade to Pivot Posi-tilt face rest base
- PTB** - Power tilt back



### Standard Upholstery Color Choices:



### Upholstery Upgrade Natural Touch Choices:



## Caring for your table:

In order to maintain your product's quality and ensure the life of your product, we suggest these preventative steps to protect the upholstery fabric on your table and accessories.

- Avoid storing your product in extreme cold temperatures.
- Excessive cold may cause the upholstery to crack.
- Guard your upholstery fabric against heat. Extreme heat may cause the upholstery to stretch and pucker. Keep away from heaters, electrical appliances, and direct sunlight.
- Protect your table and accessories from oil and makeup stains with sheets and fleeces.
- Clean your upholstery fabric every day.



## Upholstery Cleaning:

- For day to day cleaning, a solution of 10% mild household liquid dish soap with warm water, applied with soft damp cloth. Rinse with clean water and dry.
- For stains, dampen a soft white cloth in a one to one (1:1) solution of Fantastik® and water OR Formula 409® and water. Rub gently and rinse with a water dampened cloth.
- For more difficult stains, dampen a soft white cloth with a solution of household bleach (10% bleach / 90% water). Rub gently and rinse with a water dampened cloth to remove bleach concentration.
- There is a disinfecting wipe on the market called Protex Ultra Disinfectant Wipes. Many of our customers use these.
- WHAT NOT TO USE: Using the wrong cleaning agents you can void your warranty and crack, dry out and destroy your vinyl. Do not use any cleaning agents that contain alcohol, harsh chemicals or abrasives. When in doubt, give us a call!