

SOFTCAIRE DX (DELUXE SYSTEM)



ITEM #	DESCRIPTION
H88410	SoftCaire DX (Deluxe System)

Quick Connector with Cap

- Ensures Mattress remains inflated during power outage or patient transport.

CPR function

- Allows rapid deflation in case of medical emergencies.
- Designed for one handed operation.
- Positioned at head end.

Pillow Function

- Head zone comprised of three fully inflated static air cells to provide patient comfort.

Square Cell-on-Cell Bladder

- Unique two-layer cell structure provides both alternation therapy and constant support in case of power-outage or transportation.

Heel Function

- Provides single air cell deflation to float the heels with zero pressure.

Spring-loaded Tube

- Spring are emplaced within air tubes to prevent unnecessary bending of tubes and ensure consistent airflow.
- Incorporates an 8 inch cell-on-cell mattress design that can maintain basic support during power failure or patient transport.
- Utilizes Bi-elastic Nylon/PU top cover that is water resistant, vapor permeable, and anti-bacterial.

SPECIFICATIONS

Weight Setting Control:	Digital
Dimension:	12.2"L x 4.9"W x 8.25"H
Weight:	5.5 lbs
Air Output:	9 liter/min
Case Material:	Flame Retardant Abs
Pressure Range:	20-55 mmHg
Cycle Time:	10/15/20/25 min
Power Supply:	120V/60 Hz
Inflated Dimension:	80"L x 35"W x 8"H
Alternating:	A/B Alternating
Structure:	20 x 8" Cell on Cell
Cover Material:	Nylon/PU Quilted with Zipper
Cell Material:	Nylon/TPU
Base Material:	Durable Nylon/PVC 840D
Maximum Weight Capacity:	450 lbs
Flame Retardant Standards:	BS6807 Crib 5/ BS7175 Crib 5

SOFTCAIRE DX (DELUXE SYSTEM)



18 MONTH LIMITED WARRANTY

Your Rhythm product is warranted to be free of defects in materials and workmanship for eighteen (18) months of the original consumer purchaser. This warranty does not cover device failure due to misuse, negligence, improper storage or handling, improper operation and unauthorized modifications. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, which are subject to normal wear and need periodic replacement.

For warranty service, it is recommended but not required that the product be returned to the Rhythm Healthcare dealer through whom it was originally purchased. However, all Rhythm Healthcare dealers are qualified to assist you in obtaining warranty service. If the product is to be returned to Rhythm Healthcare, prior authorization will be required.

This product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have a question about your Rhythm Healthcare product or this warranty, please email us at contactus@rhythmhc.com