

Serene Elite

Alternating Pressure Redistribution System with Low Air Loss

- Suitable for patients with high to very high risk for pressure ulcers
- Ventilated low air loss for patients' upmost comfort
- Cell-on-Cell bladder to avoid bottoming out
- Unique heel relief function for comprehensive care for patients
- The design of five different modes of treatment fulfilling various patients' needs













Apex Medical Corp. Serene Elite



Serene Elite

Alternating Pressure Redistribution System with Low Air Loss


Pressure Relief During Alternating

Serene Elite

Pressure mapping reveals that in alternating mode, bladder A & B take turns to inflate and deflate can significantly relieve the body pressure

5 Modes to Fulfill Individual Needs

Max Firm
Alternating



4-Cycle Time (5/10/15/20 mins)

10 Comfort Levels



Heel Relief Function Single air-cell deflation can float the heels with zero pressure



Visual and Audible Alarm

- Low Pressure
- Power Failure
- Tech Support

High Compressor Volume Airflow





Specifications are subject to change for improvement without notice

Ordering Information (For North America Only)

Model		Serene Elite		
		36"	42"	48"
System	With Quilted Cover	PD04501	PD04504	PD04507
	With Non-Quilted Cover	PD04502	PD04505	PD04508
	With Side Bolster, Quilted Cover	PD04503	PD04506	PD04509
Mattress	With Quilted Cover	PM09407	PM09501	PM09601
	With Non-Quilted Cover	PM09408	PM09502	PM09602
	With Side Bolster, Quilted Cover	PM09409	PM09503	PM09603
Pump			PD04001	

4-Way Stretch PU / Nylon PU

- Low friction and shear
- Vapour permeable
- High-flame retardant
- Water resistant





Zipper-all-around Cover Reduce ingress of unwanted

remove for cleaning

fluid and simply one step to













Quick Connector Once with cap on, instantly switch to static air mattress for transportation need