



Links Medical Products, Inc.
9247 Research Drive
Irvine, CA 92618
Tel: (949) 753-0001
Toll Free: (888) 425-1149
www.linksmed.com

*At Links Medical, we are your link to
Innovative Caring Products*

ALL NEW SILENT KNIGHT[®] Advanced Pill Crushing Device

Another move forward from the market leader

The Silent Knight has long been a significant and highly trusted system component in medication management. Long term care facilities and hospitals in the United States and around the world use the Silent Knight everyday.

Thousands of healthcare professionals and tens of thousands of patients have benefited from the performance of this widely used device.

This brochure will introduce you to the advanced features and benefits of the New Silent Knight.



Links Medical Products Inc.[®]



Contoured handle



Onboard tray



Ribbed crusher plate

ADVANCED ENGINEERING & DESIGN

Advanced engineering and design components have been applied to the original Silent Knight's basic operating and crushing features to create the New Silent Knight.

The patented dual action "crushing and grinding" movements of the New Silent Knight result in silent, smooth, effortless and effective crushing of pills to fine powder.

The crusher handle is designed for all hand sizes to fit comfortably. User safety is enhanced with the specially designed internal handle shock absorber. And the surgical quality stainless steel crushing components, along with a solid aluminum base and handle, make the New Silent Knight a powerful, highly effective and durable piece of medical equipment.



Crusher plate action

EXCLUSIVE PILL CRUSHER POUCH

The Silent Knight® does not function alone. Its companion item the patented Silent Knight® Pill Crusher Pouch is essential to guarantee maximum performance while reducing cost, waste, malfunction and lost time.

The pouch's distinctive curved bottom structure assures full recovery of crushed medication. The thick, durable construction effectively reduces pharmaceutical cross contamination and eliminates loss of medication.

Note: No other pill crusher pouch will effectively and safely function in the Silent Knight®. Use of any other pill pouch will compromise product performance and outcome.



Step 1

Step 2

Step 3

IN-SERVICE STEPS

Step 1. Place pill(s) into Silent Knight Pill Crusher Pouch. Lift handle and place pouch down between the crusher plates to initiate crushing.

Step 2. Press handle down to crush pill(s). Lift handle and pull pouch upward to half way up the crusher plates.

Step 3. Then press the handle down to activate grinding movement. Press handle down several times to grind pill(s) into fine powder. Repeat if needed.

Step 4. Remove the pouch and add the medication to an appropriate serving medium for administration.

Visit www.linksmed.com/silentknightinservice for a full in-service video.

LIMITED LIFETIME WARRANTY

The Silent Knight is guaranteed for life to be free of materials and parts defects. If the unit fails due to parts failure it will be replaced by Links Medical at no charge to the owner.

SPECIFICATIONS & ORDERING INFORMATION

Silent Knight Pill Crusher
Materials: Stainless steel, anodized aluminum. Latex free.
Item Number SK0500

Silent Knight Pill Crusher Pouch
Materials: Low Density Polyethylene (LDPE).
FDA approved for pharmaceutical packaging.
Item Number PC1000

ACCLAIMED BENEFITS OF THE NEW SILENT KNIGHT®

- SILENT** Won't disturb patients or staff.
- ERGONOMIC HANDLE** Fits all hand sizes.
- EFFECTIVE** Powerful dual action crushes pills to powder in a few seconds.
- PATIENT SAFETY** Eliminates cross contamination.
- USER SAFETY** Maximum results with minimal effort. Contour cover encloses internal moving parts.



Protected by US patents: 7,051, 963 B1, D497, 543 S, 5, 915,637, 9,125,796.
All trademarks are the property of Links Medical Products, Inc.
© 2016 All rights reserved.