


# Gel Foam X-Cushion




Skil-Care  
Design

Order #	Description	Size	Unit
751026	Gel Foam X-Cushion - Vinyl Only	18" x 16" x 2.5"	Pair
751027	Gel Foam X-Cushion - Vinyl w/Low Shear I Cover	18" x 16" x 2.5"	Pair
751030	Gel Foam X-Cushion - Vinyl Only	16" x 16" x 2.5"	Pair
751031	Gel Foam X-Cushion - Vinyl w/Low Shear I Cover	16" x 16" x 2.5"	Pair
751032	Gel Foam X-Cushion - Vinyl Only	20" x 16" x 2.5"	Pair
751033	Gel Foam X-Cushion - Vinyl w/Low Shear I Cover	20" x 16" x 2.5"	Pair
751034	Gel Foam X-Cushion w/Coccyx Cut-Out and Low Shear I Cover	16" x 16" x 2.5"	Pair
751035	Gel Foam X-Cushion w/Coccyx Cut-Out and Low Shear I Cover	18" x 16" x 2.5"	Pair
751036	Gel Foam X-Cushion w/Coccyx Cut-Out and Low Shear I Cover	20" x 16" x 2.5"	Pair

**Available Low-Shear** cover reduces skin friction and shear. Fluid-proof. Wipes clean for easy care.

**Vinyl cover** is incontinent-proof and wipes clean for easy care.

**High-density foam** base prevents bottoming-out and increases sitting comfort.



**Straps** with buckle secure cushion to wheelchair.

**Water-based gel** surface distributes weight to relieve pressure. Its cooling action reduces skin damaging perspiration.

**New X-Design** creates four separate gel sections for maximum stability.

## Purpose

Skil-Care's Gel Foam X-Cushion is a combination gel and foam cushion designed to relieve pressure. The waterbased gel surface provides cooling action that reduces skin-damaging perspiration.

## Important Safety Reminder

The Gel Foam X-Cushion with a cloth outer cover has straps with buckles that must be secured to the wheelchair to prevent the cushion from sliding off the seat.

## Easy Care

The vinyl cover may be wiped clean using a damp cloth and spray-type cleanser. A mild disinfectant may be used. The Low-Shear outer cover can be laundered and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach.