

NEW



TAP PAP NASAL PILLOWS

No headgear necessary!



Mouthpiece molds to the upper teeth to provide a custom and secure fit even if the bottom teeth are not in contact with the mouthpiece.

Unique adjustment joint with nut & bolt connection allows the nasal pillows to be set to the desired angle

Provides a stable fit regardless of the sleep position

Patients report not being able to feel the mask on their face.

3 sizes of nasal pillows
(small, medium and large)

90 day compliance guarantee for new users of the mask, 30 day guarantee for experienced users

Order#

CM024

Mask kit with Small, Medium, and Large Pillows and optional headgear

CUCM024S Small Pillows

CUCM024M Medium Pillows

CUCM024L Large Pillows