

# Suspension Unit

The Suspension unit is a carefully designed harness that is comfortable to wear both sitting and standing and is easy to put on/take off. In addition, the unit is adjustable to provide maximum stability and support to the pelvis. The unit is washable with soap and water or standard disinfectants.



## **Support**

The back side of the unit is a flexible gel pad shaped like a bicycle seat. The front support is a wide pad of compressible foam.

## **Putting On/Taking Off**

The unit is opened/closed with large, easy to use snaps.

## **Adjustment**

The unit is adjusted for fit and comfort by 4 easy to use "tug and release" slides.

## **Connection**

Connecting the Suspension Unit to the Mobility Unit is done by lifting a ring over a ball. The simple connection can be done with a care giver or by the user. (Depending on the user's capability)