



Wrist & Thumb Positioning System

The **Wrist & Thumb Positioning System** is the most versatile and flexible system of its kind. Prescribe it as a stand-alone functional alignment aide, or use it to extend a TheraTogs strapping application from the shoulder or elbow to the hand.

No other wrist or hand support product offers the adaptability, flexibility, and patient comfort of a TheraTogs system!

The hook-sensitive wrist and thumb unit comes in Pediatric (S, M-L) and Adult (S, M, L) sizes for either the left or right hand. Made of the same foam-backed material as our TheraTogs garments, the Wrist & Thumb Positioning System is both comfortable and effective in supporting and maintaining joint position gains that are available with manual correction.



Indications/Applications:

Use the Wrist & Thumb Positioning System to:

- q Improve wrist and thumb stability
- q Restore the palmar arch
- q Reduce flexible ulnar deviation
- q Assist in reducing forearm pronation or supination
- q Relieve pain associated with malalignment




Wrist & Thumb System Contents:

The Wrist & Thumb Positioning System comes with these components:



Improve wrist extension and radial deviation

- q One Wrist & Thumb System (*specify Left or Right when ordering*)
- q One Forearm Supination Wrap
- q One Rotation Strap
- q One Limb Cuff
- q Three Split Straps
- q Size-appropriate tabs
- q Marker dots
- q Printed, illustrated instructions for donning, fitting, and strapping.



Help to restore the palmar arch.

If you can correct wrist and thumb alignment with your hands, without force, you can use TheraTogs to carry over that correction – positioning your client's day-long activities in corrected functional alignment.