

AUTO WAVE 695

Radial Pressure Pulse Therapy



Mettler Electronics Corp.

AUTO WAVE 695

Radial Pressure Pulse Therapy



Features

- Silent built-in compressor
- 7-color touch screen
- Ergonomic hand piece
- Treatment protocols with images
- Multiple language capability



Indications

- Helps relieve minor muscle aches and pains
- Temporarily increases local blood circulation



Applications for:

- Orthopedics
- Physical Therapy
- Podiatry
- Chiropractic





Radial Pressure Pulse Therapy is also known as Extracorporeal Pulse Activation Therapy. The Auto*Wave 695 has a compressed air ballistic pulse energy generator that transfers pulse energy to the transmitter in the hand piece.

The motion and weight of the projectile, which is accelerated by compressed air, produces kinetic energy that is converted into pulse energy when the projectile strikes the surface of the transmitter.

The acoustic pulse energy is transformed into impact energy from the contact area into the therapeutic target tissue radially.

☼ Specifications

- **Pressure Range:** 1.0 - 4.0 Bar (0.1 Bar steps)
- **Frequency:** 1 - 17 Hz (1 Hz steps)
- **Pre-Set Protocols:** Yes
- **User-Defined Protocols:** 5,000
- **Dimensions:** 15.5" (L) x 15.3" (W) x 7.6" (H)
- **Weight:** 18.7 lbs (Headpiece excluded)



Ballistic Mechanism of Radial Pressure Pulse Therapy

☼ Transmitter

Transmitters						
Type	6 A (D = 6mm Acupuncture)*	15 R (D = 15mm Radial)	15 F (D = 15mm Focus)*	15 D (D = 15mm Deep)*	20 R (D = 20mm Radial)	36 R (D = 36mm Radial)
Penetration Depth	0 - 40 mm	0 - 40 mm	0 - 20 mm	0 - 60 mm	0 - 50 mm	0 - 50 mm

* Optional Accessories

Certification

CE 0537

CE (MDD)



Health Canada



New Custom Cart

- Four locking wheels
- Large handle to easily move the cart
- Three internal shelves to store accessories and supplies
- Easy opening magnetically secured door



Mettler Electronics Corp.