

Alliance™ Patient Lifts

Alliance™ Patient Lifts, Slings and accessories provide the utmost in safety, comfort and affordability to caregivers and patients. Their ergonomic design makes transfers as simple as possible for patients who are able to bear partial weight. Alliance HE lifts feature robust construction and are easily upgradeable from hydraulic to battery-powered at an affordable price and with no tools required.

Alliance Lift Features

- Six point multi-tilt spreader bar for additional stability and comfort
- Only lift manufactured with lift/lower emergency backup system (REF: 1902 Battery Powered Lift only)
- Interchangeable battery systems feature audible low battery warning
- Foot pedal base opening for wheelchair/commode access
- Large, easy-to-grip handles for superb maneuverability
- Low 4.5" (11 cm) base height provides maximum stability and ensures fit under lower-clearance beds
- 400 lb (181 kg) safe working load
- Also available for bariatric patients - each with 600 lb (272 kg) safe working load
 - » Heavy Duty Lift (REF: 1904)
 - » Stand-Assist Heavy Duty Patient Lift (REF: 1907)



Alliance HE Patient Lift
HCPCS E0630†
 1900 – Hydraulic Patient Lift HE



Alliance Battery Powered Patient Lift
HCPCS E0635†
 1902 – Battery Powered Patient Lift

Alliance HE Lift quickly converts to battery power!



HCPCS E0635†
 1908 – Battery Powered Patient Lift HE
 1910 – Battery Power HE Upgrade Kit
 (for use with REF: 1900)
 19526 – Replacement Battery



Alliance Heavy Duty Patient Lift
HCPCS E0635†
 1904 – Heavy Duty Battery Powered Patient Lift (Bariatric)

Alliance Digital Scale

Includes all necessary hardware for easy installation. Maximum 600 lb (272 kg) calibrated weight capacity. Compatible with REF: 1902, and 1904.



19550 – Alliance Digital Scale



Alliance Stand-Assist Patient Lift
HCPCS E0635†
 1906 – Stand-Assist Patient Lift
 1907 – Stand-Assist Patient Lift HD (not shown)
 600 lb (272 kg) safe working load

Patient Transfer Sheets (not shown)

Low friction nylon fabric allows for easy patient repositioning and is specially treated to guard against static buildup.

- 19589 – Small
28.5" x 27.5" (72 cm x 70 cm)
- 19590 – Medium
42.5" x 31.5" (108 cm x 80 cm)
- 19591 – Large
53" x 39.5" (135 cm x 100 cm)