


HOW MUCH WILL BE COVERED?  
= GREEN AREA




HOW CAN I MOVE THE LIFT?  
= FROM A TO B



HOW CAN I MOVE THE LIFT?  
= ROOM COVERING SYSTEM



HOW CAN I MOVE THE LIFT?  
= FROM ROOM TO ROOM (R-T-R)



# Bedroom Ceiling Mounted

- TL8-03028 Bestcare Luna 275kg / Bayonet (1)
- TL8-03101 Pendant/1320mm/ bottom exit (1)
- TL3-00702 Charger/holder / Luna handset (1)
- TL8-01968 BestLift Luna Sbar/strap/275kg (1)
- TL6-01022 Track 75mm/4M (1)
- TL6-01132 End Stop 75mm Track (2)
- TL6-01192 End Covers 75mm Track (1)
- TL6-02002 Cart for 75MM Track Bayonet

