


RED LIGHT THERAPY DEVICES

PROFESSIONAL MANUFACTURER IN CHINA SINCE 2015



WHY NEED RED LIGHT THERAPY?

Red light therapy, also known as low-level laser therapy (LLLT) or photobiomodulation, offers a range of benefits that make it a popular choice for various health and wellness concerns. This non-invasive treatment utilizes specific wavelengths of red light to penetrate the skin and stimulate cellular repair and rejuvenation processes.

It can help reduce inflammation, improve collagen production, and enhance overall skin tone and texture, making it effective for addressing issues like acne, wrinkles, and scars. Furthermore, red light therapy is known for its pain-relieving properties and its ability to accelerate wound healing. By increasing circulation and reducing inflammation, it can aid in managing chronic pain conditions and promoting faster recovery from injuries.

Additionally, red light therapy has been shown to boost mood and improve sleep quality by regulating circadian rhythms and promoting the production of serotonin, often referred to as the "feel-good" hormone. Whether you are looking to enhance your skin's appearance, alleviate pain, or improve your overall well-being, red light therapy provides a safe and effective solution with a wide range of applications for various health concerns.



Welcome to BioMol.

We are one of the top-rated manufacturers of Red Light Therapy devices, with a 2200 square meters factory located in Shenzhen China since 2015.

With over 60 staff members, we ensure a production capacity of 2000pcs per month for quick lead time. Our products are of top quality and built with CE, RoHs, FCC and FDA certificates. Furthermore, our factory complies with ISO 9001 standards.

For an enterprise, innovation is the foundation of long-term development. BioMol Therapy has a strong R&D team of over 10 people to ensure continuous product innovation each year. All our products are patented with unique designs. We take pride in being the ONLY company producing Red Light Therapy devices with Voice Control. Additionally, we are the ONLY company to divide our product line into different series based on their applications.

Our product series include:

Portable Series

Household Series

Household Pro Series

General Series

Beauty Medical Series

Fitness Series

RP Series

RL Series

Wearing Series

Pets Care Series

Our products are with the following unique features:

Smart Voice control

2-8 channels with dimming & timing & pulse

Pre-set modes for quick use & customized modes with settings

The highest irradiance devices available on the market

With our professional and strong R&D team, you are always welcome to customize your exclusive products through our OEM&ODM service, and we provide FREE design service.

Choosing BioMol as your business partner will provide you with high-quality guarantee services for your market competition and long-term development.

Let's begin a win-win cooperation now!

CFO

Ivy / Miss Chao

WHY CHOOSE US?



Professional manufacturer with 16 years experience in red light therapy indutry.



Strong R&D team with 20+ members with 3-4 new series annually.



Unique design with voice control function (the only supplier on the market).



Free OEM&ODM design service.



High quality of 3-year warranty with new replacement.



FDA, CE, RoHS certificates for guarantee usage.


OUR COMPANY













CERTIFICATION



BIOMOL®

PORTABLE SERIES BM-SLIM



Skin care & anti-aging

Promote collagen synthesis: red light therapy can promote the synthesis of collagen in the skin, making it more firm and elastic.

Improve pigmentation: red light therapy can regulate collagen production, reduce pigment deposition, and improve pigmentation problems.

Shrink pores: red light therapy can promote the metabolism of skin cells, clean and shrink pores.

Just enjoy it!



6 Default Modes



6 Spectra



120 x 0.5W LEDs



Aluminum Housing



Touch Key



Small Size & Portable



BioMol Slim













Skin rejuvenation

- Spectrum: Red 630nm&660nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance





Skin whitening

- Spectrum: Red 630nm&660nm Yellow 590nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance





Acne treatment

- Spectrum: Blue 450nm /Red 630nm& . 660nm/NIR 850nm&940nm
- Timing: 15 minutes
- (5mins Blue + 10mins Blue/Red/NIR)
- Usage: once a day, 10-30cm distance





Skin firming

- Spectrum: Red 630nm&660nm NIR 850nm&940nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance





Freckle treatment

- Spectrum: Yellow 590nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance







- Oil control & purifying Spectrum: Blue 450nm/Yellow 590nm/Red 630nm&660nm/ NIR 850nm&940nm
- Timing: 15 minutes
- (5 mins B/Y/R/NIR+ 10 mins Y/R/NIR)
- Usage: once a day 10-30cm distance

TOUCH KEYS


Simple to control the device






BM-Slim





IRRADIANCE TEST



STAND



Table stand



Adjustable floor stand (Optional)

PACKAGE CONTENTS

- BioMol LED red light therapy panel x 1
- Readable glasses x 1
- DC adapter x 1
- BioMol user manual x 1
- Table stand x 1









