

COLD RUSH®

INDICATIONS FOR USE

- For cold therapy pain relief, as needed

FEATURES AND BENEFITS

- Consistent cooling and excellent compression, even at high elevation.
- Brushless motor produces a whisper quiet unit that allows patients to get the rest they need to properly rehabilitate.
- Highest ice-to-water ratio in the industry and a design that produces less heat transfer from the motor allows therapy pain relief for up to 6 continuous hours.
- Unique locking lid prevents spills
- Insulated walls improve cooling



COLD RUSH® COMPACT

INDICATIONS FOR USE

- For anything that may benefit from cold therapy

FEATURES AND BENEFITS

- Convenient, quiet and easy-to-use
- Single flow setting
- Easy fill instructions



ORDERING INFORMATION

| Description | Cold Rush Part # | Cold Rush Compact Part # |
|------------------------------|------------------|--------------------------|
| Device Only USA | B-232000010 | B-232100010 |
| Device w/ Knee Pad | B-232000011 | B-232100011 |
| Device w/ Unv Pad | B-232000012 | B-232100012 |
| Device w/ Hip Pad LT | B-232000013 | B-232100013 |
| Device w/ Hip Pad RT | B-232000014 | B-232100014 |
| Device w/ Shoulder Pad SM | B-232000015 | B-232100015 |
| Device w/ Shoulder Pad LG | B-232000016 | B-232100016 |
| Device w/ Ankle Pad | B-232000017 | B-232100017 |
| Device w/ Rectangular Pad SM | B-232000018 | B-232100018 |
| Device w/ Lumbar Pad | B-232000019 | B-232100019 |

ACCESSORY

| Description | Part # |
|------------------------|-------------|
| Replacement Power Cord | B-301090001 |