

Comfy™ Deviation Finger Extender

For the treatment and prevention of ulnar and radial deviations

DEVIATE

Pivot/swivel mechanism adjusts using hex wrench (included) to accommodate for ulnar and radial deviation

ADJUSTABLE

Patented reinforced spine bends for individualized wrist and forearm needs

PROGRESSIVE FINGER EXTENSION

Two graduated and removable rolls gradually increase extension of the MCP and PIP joints

COMFORT

Super soft, antimicrobial, removable Headliner or Terrycloth cover is machine washable



High Reimbursing
PDAC code
L3915 or L3916

DF-101 (h)	Deviation Finger Extender (h= headliner)
DF-101-AS (h)	Deviation Finger Extender in Adult Small (h=headliner)
PDF-101- (S,M,L)	Pediatric Deviation Finger Extender

** Should be used only in connection with a care plan by a trained health care professional.