


PERFECT FOR PATIENTS
PERFECT FOR CAREGIVERS

Enhance

Enhance Mattress System is a low to medium risk system for patients who are at risk of developing Pressure Ulcers according to your sound clinical judgment. Enhance can also be used for patients who have an existing stage 1, 2 and 3 Pressure Ulcer, in conjunction with your hospital's policy on pressure area management.



Master Control Unit

- Dimension: (L)10.1" x (W)4.5" x (H)8.3"
- Weight: 5.1 lbs
- Rating Voltage: 110-120 volts
- Max Current: 0.1A max.
- Rating Frequency: 60Hz
- Cycle Time: 10 mins
- 2 in 1 alternation
- Min/Max Pressure: 16-32 mmHg +/- 6 mmHg
- Flow-rate: ≥ 4.5 l/min
- Protection Type: Class I, Type BF
- Low Pressure Alarm: Visual / Audio

Mattress Overlay

- Model No: Enhance
- Dimension: (L)80" x (W)36" x (H)5"
- Weight: 16.7 lbs
- Weight Capacity: 350 lbs
- Cell Material: Blue Nylon with PU Coating
- Cell Number: 20 cells
- Cover Material: Blue Nylon with PU coating finish
- Base Material: Black Polyester with PVC Coating
- CPR Quick Release System

CREATING INNOVATION AND DIFFERENCE



Enhance




Enhance

provides you, the caregiver, with peace of mind and your patient, with world class care & comfort

Toll Free:
877-367-5507

LA, CA OFFICE
TEL: 626-275-4063
TAMPA, FL OFFICE
TEL: 813-854-5464

Perfect for patients and caregivers



PRESSURE MAPPING

Enhance provides constant low pressure by redistributing the patient's pressure over a greater area than a standard mattress.

• The Enhance overlay mattress can accommodate patient weights up to 350 lbs

• Interchangeable air cell design provides easy cleaning and minimizes the proliferation of bacteria

• The Enhance, with it's 5 inch high cells, provides safe support and personal comfort to patients

• Modular design ensures ease of use and cleaning



Enhance is your partner in protection and treatment of pressure ulcer




Cable management keeps electrical cables off the floor



CPR quick release system deflates the system within 20 seconds

Low slip base cover is a plus to keep the mattress in place on the bed



Effective micro low air loss for moisture management



Slimline pump design is unobstructive

Adjustable control allows for 10 levels of comfort depending on the patient's therapy requirement.



The Low-Pressure Alarm is triggered when the air pressure of the mattress is lower than the default setup.

2-in-1 alternating air cycle ensures the maximum protection.



Hanging bracket fits on all type of beds