

The Dura-Gel™ SP III (Skin Protection) wheelchair cushion reduces surface pressure through the use of non-aqueous Dura-Gel™ placed in three pressure zones. The cushion uses a high-density base foam in combination with a visco-elastic topper for added comfort and also positions the user with a moderate slope from front to back.



- ☐ For use with moderate risk patients
- ☐ Reduces tissue interface pressures at critical points
- ☐ Features a specialized Dura-Gel[™] center layer surrounded by comfortable foam layers on the top and bottom of the cushion
- ☐ The soft, stretchable lightweight cover reduces shear/friction and is fluid resistant and has a strap at the back to secure the cushion firmly to the chair
- ☐ Accommodates a wide range of body shapes, sizes and weights
- □ 350 lb maximum weight capacity
- ☐ Lightweight
- ☐ 18 month limited warranty



Designed to meet the requirements of HCPCS Code: E2603

8930166 16" x 16" x 3" 8930186 18" x 16" x 3" 8930188 18" x 18" x 3" 8930206 20" x 16" x 3" 8930208 20" x 18" x 3"

Designed to meet the requirements of HCPCS Code: E2604

8930228 22" x 18" x 3" **8930220** 22" x 20" x 3" **8930248** 24" x 18" x 3



