

Alternating Pressure Systems

Genesis I and II Series Alternating Pump and Pad System

The Genesis Series Pressure Pump and Pad System is designed with advanced technology that alternately inflates and deflates air cells in the pressure pad. The small compact pump conveniently attaches to the bed and is available with fixed or adjustable pressure. The Genesis Series features an output of eight liters per minute with a five minute cycle time and a one year pump warranty.



APP-1000-GBN



Genesis I Pump



Product Features

- Small compact design with reinforced hooks located on back of the pump
- Pump produces 8 LPM airflow with five minute cycle time
- Bubble cell design delivers comfort while providing superior therapeutic treatment
- Pad is made of long lasting durable heavy gauge premium vinyl
- Variable pressure setting for optimum comfort and support (APP-200-GBN model only)
- Quiet pump alternately inflates and deflates the air cells
- 90 day warranty on mattress pad
- One year warranty on pump
- HCPCS Code: E0181* (pump), A4640* (replacement pad)

Product Specifications

- Power: 120 VAC
- Input Frequency: 60 Hz
- Pressure Range: 80 mmHg (APP-1000-GBN), 40-90 mmHg (APP-2000-GBN)
- Air Volume: 8 liters per minute
- Cycle Time: 5 min/60 Hz
- Pump Weight: 3.3 lb.
- Pad Weight: 4 lb. (APP-BPD), 4.8 lb. (APP-DLXP)
- Pump Dimensions: 10.2" x 5" x 3.1"
- Pad Dimensions: 78" x 35", with tubing
- Cell Depth: 2.5"
- Power Cord: 10 ft.
- Weight Capacity: 250 lb. on bubble pad



APP-2000-GBN

Whenever bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce the risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use.

Item #	Description	Pad Type	Qty
APP-1000-GBN	Genesis I Fixed Pressure Pump and Pad System	Basic	each
APP-2000-GBN	Genesis II Alternating Pressure Pump and Pad System	Deluxe with end flaps	each