




EZ Way Hourglass Slings

EZ Way, Inc.


"Your Total Patient Lift Solution"

The hourglass-shaped sling is designed for maximum comfort including special head support. Available with solid or mesh material and 'toileting hole' option.




- 500235 - Small
- 500103 - Medium
- 500211 - Large
- 500187 - XL
- 500172 - XXL

Hourglass Sling



Comfortable, Strong & Durable




- 500267 - Small
- 500268 - Medium
- 500269 - Large
- 500270 - XL
- 500259 - XXL
- 500258 - XXXL

Hourglass Sling w/ Toileting Hole


Common Applications:

- Used with your Patient/Resident where the sling would be left underneath them while sitting in a chair/wheel chair.*
- Amputee above the knee, either single or double.
- Patients with sensitive skin, brittle bones and other special conditions which might apply.
- Patient/Resident wearing a dress or skirt and dignity prevails.
- Not for use with combative or uncooperative Patient/Resident.




- 500126 - Small
- 500061 - Medium
- 500168 - Large
- 500175 - XL
- 500203 - XXL

Mesh Hourglass Sling



Breathable, Strong & Comfortable



- 500254 - Small
- 500251 - Medium
- 500252 - Large
- 500253 - XL

Mesh Hourglass Sling w/ Toileting Hole

*The hourglass sling is manufactured in different materials with special options. If this is to be left under a Patient/Resident for a period of time, EZ Way recommends using the mesh hourglass sling. Ask your EZ Way Territory Sales Manager or call 800-627-8940 for more information.