# **RAV IV Laser Therapy**

# **INFINIUM**®



Laser Focused Healing at the Speed of Light CLINICAL STRENGTH THERAPY SYSTEM

# LASER THERAPY

# Infinium RAV IV Laser Deep Penetrating Therapy

The Infinium RAV IV is a clinical strength Class 4 laser therapy system that induces cellular level reactions to promote the healing and recover of soft tissue. The RAV IV emits powerful coherent laser in specific wavelengths. This deeply penetrates beyond the skin to underlying tissues, including muscles, tendons, and ligaments. The RAV IV reduces inflammation by influencing the release of pro-inflammatory cytokines and promoting vasodilation (expansion of blood vessels), which enhances blood flow and reduces swelling. The RAV IV also stimulates the production of mitochondrial activity, increasing ATP (adenosine triphosphate) production. ATP is crucial for cellular repair and regeneration processes.

- >> 810nm / 980 nm Wavelengths
- >> Up to 15 Watts Energy Output
- Continuous Wave or Pulsed Modes
- Easy to Operate Touch Screen
- Preprogramed Therapy Protocols



- Multiple Hand-Pieces
- >> Painless Non-Invasive Treatment
- Low Cost of Ownership



# Infinium RAV IV Therapy Applications

#### **Head and Neck**

- Cervical Disc Herniation
- · Cervical Radiculopathy
- · Cervical Soft Tissue Strain
- · Cervical Spondylosis
- Cervical Stenosis
- · Headaches, Cervicogenic
- · Headaches, Migraine
- · Headaches, Tension Sinusitis
- TMJ Syndrome
- · Trigeminal Neuralgi

#### **Back & Shoulder**

- Back Lower, Upper, Lumbar
- · Back Disc Pathologies
- · Lumbar Facet Joints
- · Lumbar Spondylosis & Stenosis
- · Lumbar Radiculopathy
- · Sacroiliac Joint
- Tailbone
- Sciatica
- · Rotator Cuff
- · Shoulder Frozen, Strain, Sprain

#### **Upper Extremities**

- · Biceps, Tendonitis
- Carpal Tunnel Syndrome
- · Epicondylitis Lateral, Medial
- Fractures
- · Joint, Wrist
- Olecranon Bursitis
- · Triceps, Tendonitis
- Ulnar Nerve Pathologies
- · Joints, Fingers
- Soft Tissue Injury

#### Hip & Legs

- · Abductor Ilio-Psoas / Band
- · Bone Leg Fracture
- Hamstrings
- · Hip Joint, DJD, RA, Periarticular
- · Knee Sprain, Bursitis, DJD, RA
- Ligament Cruciate Ant/Pos
- Ligament Meniscus Injury/Sprain
- · Tendonitis Infra-patellar
- Tendonitis Pat

#### **Foot & Ankle**

- Ankle Sprain I III
- Bone Facture, Stress
- Hallux Valgus
- · Heel Pain
- Plantar Fasciitis, Strain
- Shin Splints
- Metatarsalgia
- Morton's Neuroma
- Turf Toe
- Peripheral Burning, Numb
- Tendonitis Achilles, Tibial
- Tarsal Tunnel Syndrome

#### General

- · Arthritis, Neck
- · Arthritis, Shoulder
- · Arthritis, Back
- · Arthritis, Hip
- · Arthritis, Knee
- · Arthritis, Ankle
- · Arthritis, Wrist
- · Arthritis, Fingers
- · Arthritis, Elbow
- Fibromyalgia · Wound Pain

#### **INCLUDED HANDPIECES**

#### THERAPY READY PROTOCOLS









### **Technical Specifications**

Wavelength 810nm,980nm

**Maximum Power** 15W

Operation Mode Continuous Wave, Single Or

Repeat Pulse

Pulse Duration 10µs-3s

Repetition Rate 0.2Hz-20KHz

Control Method 8 Inch Touch Screen

**Lifetime** 10,000 Hours

Warranty 1 year / service to 5yrs

**Current Rating** 110/220VAC, 5A, 50/60 Hz

Weight 4KG / 8.8lbs

**Dimensions** 31(H)\*22(W)\*26(L) cm

12.2(H)\*8.67(W)\*10.2(L) in

# **INFINIUM®**

## **RAV IV Laser Therapy Machine**

## **Flexible Lens Options**

Unique lenses allow the laser to be adapted for various treatments. Each lens focuses or disperses the laser light for specific conditions or areas of the body to optimize treatment.

The included kit provides contact and non-contact lenses from 5 mm to 30 mm. These lens options adjust the depth of penetration and the intensity of the laser light, which can enhance patient safety, maximize therapy results, and minimize potential side effects.





