



AQUATEC

Aquatec[®] **Marine**

Bathing simply reinvented



The *Aquatec Marine* makes bathing simple and has been designed with **comfort, ease of use** and **safety** in mind. Weighing only 22 lbs., the *Marine* is **lightweight** yet able to support up to 308 lbs. Having been tested against the toughest safety standards, the bathlift has a robust and stable structure allowing you to bathe with peace of mind.

- Easy to use with a smooth riser mechanism
- Easy to assemble and maintain
- Lightweight and easy to move and handle



Features and Options



Floating handset

This easy to use, floating hand control has been **ergonomically** designed and can be placed in its handset holder when not in use.



Open base plate

There is no need to remove the *Marine* from the tub when you clean thanks to the open base plate.

Comfortable bathing experience

While straightforward in nature, the *Marine* has many features including its easy to use hand-control, smooth moving riser mechanism, its ability to lower to within 2.3" of the floor of the tub and its shower seat, allows more room to bathe.

Assembly and maintenance

The *Marine's* simple design makes assembly and maintenance very easy. Cleaning is effortless due to its smooth surface and open frame base plate, allowing you to clean the *Marine* while it is still in the bath. Its lightweight design makes it easy to move in and out of the bath, all of which allow for a **easy, comfortable** and **safe** bathing experience.

Technical data

For more information about this product, including the product's user manual, please see our website, www.clarkehealthcare.com.

Photos showing optional white covers.

	<i>Marine</i>
Seat width with side flaps folded in	14.4"
Seat width with side flaps folded out	27.8"
Seat depth	17.5"
Seat height at highest position	16.3"
Seat height at lowest position	2.3"
Backrest width	14.4"
Height of backrest	25.7"
Base area	22.2" x 10.3"
Total weight without battery	22 lb
Max. load	308 lb.
Backrest can be reclined by up to	40°
Battery	Lead storage battery



Ph: 888-347-4537
 Fax: 734-695-2922
info@clarkehealthcare.com
www.clarkehealthcare.com

Clarke Health Care Products
 7830 Steubenville Pike
 Oakdale PA 15071