

COLD RUSH®

INDICATIONS FOR USE

- For cold therapy pain relief, as needed

FEATURES AND BENEFITS

- Consistent cooling and excellent compression, even at high elevation.
- Brushless motor produces a whisper quiet unit that allows patients to get the rest they need to properly rehabilitate.
- Highest ice-to-water ratio in the industry and a design that produces less heat transfer from the motor allows therapy pain relief for up to 6 continuous hours.
- Unique locking lid prevents spills
- Insulated walls improve cooling



COLD RUSH® COMPACT

INDICATIONS FOR USE

- For anything that may benefit from cold therapy

FEATURES AND BENEFITS

- Convenient, quite and easy-to-use
- Single flow setting
- Easy fill instructions



ORDERING INFORMATION

Description	Cold Rush Part #	Cold Rush Compact Part #
Device Only USA	B-232000010	B-232100010
Device w/ Knee Pad	B-232000011	B-232100011
Device w/ Unv Pad	B-232000012	B-232100012
Device w/ Hip Pad LT	B-232000013	B-232100013
Device w/ Hip Pad RT	B-232000014	B-232100014
Device w/ Shoulder Pad SM	B-232000015	B-232100015
Device w/ Shoulder Pad LG	B-232000016	B-232100016
Device w/ Ankle Pad	B-232000017	B-232100017
Device w/ Rectangular Pad SM	B-232000018	B-232100018
Device w/ Lumbar Pad	B-232000019	B-232100019

ACCESSORY

Description	Part #
Replacement Power Cord	B-301090001