

All Body Active + Passive Trainer



Active Mode

15 Levels of Resistance

Passive ModeMotorized Assistance up to 60 RPM

Omni-Sense System

Automatically Changes from

Passive to Active Mode

BioFeedback

Displays Left vs Right Strength



Display Measurements

Left vs Right Symmetry

Percent Active vs Passive

Workout Summary Report

BioFeedback Graphs

User Accesibility

Spasm Detection

Wheelchair Hooks

Voice Guided Workout

Omni-Grip Handles



