INVACARE® RELIANT™ 350

STAND-UP LIFT





The Invacare® Reliant™ Series Stand-Up Lifts are an ideal product for everyday resident handling, ensuring caregiver security and resident peace of mind during the transfer. Intended for use with weight-bearing residents as well as those needing rehabilitation support, the Invacare Reliant Stand-Up lift will be a valuable addition to your safety management system.

Features

Innovative adjustment features adapt to a wide variety of body sizes and shapes

Low friction casters significantly improve rollability and contribute to caregiver safety

Multi-function slings allow for quick toileting, stand assist or full-support seated transfer

Adjustable leg support features a soft pad with lateral contour that secures legs in three positions for comfort at any height

Wide, non-slip footplate provides stable base for transfer

Under-bed clearance of 4.5" allows accessibility to most nursing home, hospital and homecare bed styles

Available with power opening base (Model # RPS350-2)

Removable battery eliminates downtime

Ergonomic features make it ideal for toilet transfers, as well as for dependent and rehabilitation purposes

Features include fast, simple sling attachment that protects against accidental disengagement; wide push handle; good access for clothing removal; and comfortable slings with extended back support that prevent pressure in the under-arm

Adjustable leg strap option secures leg against knee pad for optimum positioning on footplate

Manufactured in an ISO13485 facility







- 1. Electronics
- 2. Extra-large standing sling* Model # R131
- 3. Extra-large transport sling* Model # R136

*Complete instructions for lift usage and sling procedures are found in the Invacare operating manual. Only Invacare® slings and accessories should be used on Invacare® lifts.



USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 668-2337

Invacare Corporation www.invacare-cc.com

Canada

994 Hargrieve Road London, Ontario N6E 1P5 Canada (800) 668-2337

©2013

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 13-178 130306















| Model # | Height Base At Sling Hook Up Dimensions | Casters | Weight | Weight Capacity | Low Battery Alarm | Emergency Controls |
|----------|---|---------------------|---------|-----------------|----------------------|-----------------------|
| RPS350-1 | 39.6"-63.7" H 4.5"H Clearance 36.8"W Open 25.7"W Closed 35.5"L | 3" Front 5" Rear | 115 lb. | 350 lb. | Audio | Stop and Down |
| RPS350-2 | 39.6"-63.7" H 4.5"H Clearance 35.8"W Open 25.7"W Closed 35.5"L | 3" Front 5" Rear | 118 lb. | 350 lb. | Audio | Stop and Down |











| | Battery | Charge Time | Charger Input | Charger Output | Limited Warranty |
|----------|--------------|----------------|------------------|-------------------|---------------------|
| | | | | | |
| RPS350-1 | 2-24V DC | 6 Hours | 100-240V AC | 29.5 DC | 3 Years Lift |
| | Rechargeable | | | | 1 Year Motor |
| RPS350-2 | 2-24V DC | 6 Hours | 100-240V AC | 29.5 DC | 3 Years Lift |
| | Rechargeable | | | | 1 Year Motor |

