




Measurement Key

Length =

Circumference =

- Chest measurements are taken with patient lying down
- Pant measurements are taken with patient standing
- If in-between sizes refer to length measurements



SPIO® Upper Body Orthosis

Chest circumference is measured at nipple line.

Length is from Acromion to Greater Trochanter.

Short #/Long #

SS41/ S41	Chest Circumference: 41 cm, 23.00 cm long
SS44/ S44	Chest Circumference: 44 cm, 27.94 cm long
SS47/ S47	Chest Circumference: 47 cm, 31.11 cm long
SS51/ S51	Chest Circumference: 51 cm, 35.56 cm long
SS54/ S54	Chest Circumference: 54 cm, 39.00 cm long
SS57 / S57	Chest Circumference: 57 cm, 45.72 cm long
SS60/ S60	Chest Circumference: 60 cm, 47.00 cm long
SS64/ S64	Chest Circumference: 64 cm, 50.50 cm long
SS67/ S67	Chest Circumference: 67 cm, 56.40 cm long
SS70/ S70	Chest Circumference: 70 cm, 60.50 cm long
SS73/ S73	Chest Circumference: 73 cm, 64.30 cm long
SS77/ S77	Chest Circumference: 77cm, 68.50 cm long