

## How does the technology work?

Using a combination of dynamic air cells, a non-shearing layer and 4", 4 lb memory foam allows for dynamic weight distribution, an optimal healing environment and patients feeling **weightless**. Our technology is designed for treating AND healing stage 1-4 pressure injuries and wound care with statistically-proven and unrivaled healing results.

## **The United Mattress Products**


The United Mattress will accompany a traditional size hospital bed for the patient's care. There is a small indention near the middle section to allow for the head of the mattress and bed to be raised effortlessly. There are 3 valves for the head, middle and foot section to adjust the pressure of the mattress for the comfort of the patient.

36" wide, 80" long, 9" high 500 lb weight capacity Available in XL 84"

#### Product Information

ltem	Description	Max Weight
UM1001	Std 36" x 80"	500 lbs
UM2001	Bariatric 42" x 80"	750 lbs
UM3001	Scoop 36" x 80"	500 lbs
UM4001	48" Bariatric 48"x80"	750 lbs

NOTE: All items can be made extra long at 84"





#### Possible Amputation

Nursing home facility in Madisonville: A patient had been in the facility for a few weeks with an existing wound he'd had for years, which became even worse and was potential cause for amputation. The facility used a United Mattress as a last resort and within 10 days the Dr. changed his mind and the patient has continued to heal.



- Combination of 4lb 4" memory foam and air cells allow for patient's weight to be equally distributed
- No electricity or electric pump needed
- Stays below 32 mm of Hg
- CLINICALLY PROVEN to heal stage 1-4 wounds (study available)
- Able to control the head, middle and foot section independently of one another
- · NO NEED for wedges or heel boots
- Sturdy perimeter on all models to allow for easy patient transfer and movement.

lay float heal

# The United Air 2



There are physical elements of a mattress that can offer a comforting and soothing experience for a patient experiencing for example, extended home care, hospice and the end of life process, or a patient at a long-term care facility. This mattress should be supportive and comfortable, first and foremost, be easy to use and clean by healthcare professionals and not only heal stage 1-4 wounds, but PREVENT them from forming at all.


#### Product Information

ltem	Description	Max Weight
UM1002	Std 36" x 80"	500 lbs
UM2002	Bariatric 42" x 80"	750 lbs
UM3002	Scoop 36" x 80"	500 lbs
UM4002	48" Bariatric 48"x80"	750 lbs

NOTE: All items can be made extra long at 84"

The Air 2 mattress will accompany a traditional size hospital bed for the patient's care. There is a small indention near the middle section to allow for the head of the mattress and bed to be raised effortlessly. There are 3 valves for the head, middle and foot section to adjust the pressure of the mattress for the comfort of the patient.

36" wide, 80" long, 9" high 500 lb weight capacity Available in XL 84"



### Floating Technology

The United Mattress system distributes a patient's body weight to create an **optimal** healing environment.

Support rails allow for easy patient transfer without the ballooning effect

Valves and Hosing are labeled Head, Middle and Foot and allow the user maximum comfort and adjustability Top material is non-permeable and keeps the patient cool

Middle material is made of polyester and creates a non-shearing layer

Roho Prodigy air cells are nonpowered, adjustable and reactive to movement

2 lb, 4" memory foam works with air cells for increased

pressure distribution



Bottom material is non-permeable and grips to bed frame so there is no need to strap the mattress to the frame