

WHIRLPOOLS THERMALATORS® cleatPRO® ACCESSORIES THE PERFECT SIZE COLOR HOT



CHOICES

Customize your team logo and color.

Whitehall's optional Enviro-Glaze® powder-coat finish can be a functionally colorful alternative to a standard stainless steel finish. Enviro-Glaze® color finishes are electrostatically spray-applied and cured in an excess of 400°F, bonding to the metal for a crack-, chip- and bubble-resistant surface. The unit maintains corrosion resistance and is easy to clean, providing a colorful compliment while maintaining longevity & durability.





















SPORTS MEDICINE - WHIRLPOOLS

Whirlpools help heal sports injuries, prevent muscle spasms and speed recovery time.



MOBILE

- Time-saving—no special plumbing connection
- Reversing turbines (standard) are used to pump used water out of the tank and into an available sink or drain
- Can be easily moved, facilitating space-saving convenience
- Easily transported on heavy-duty locking casters




Show your team pride with custom colors and

logos!



- Saves time when draining water
- Intended for a fixed location
- Offers direct connections to drain and water supply. Comes standard with high capacity gravity drain





SPORTS THERAPY PRODUCTS



Mobile and Stationary Units available in: 4 packs, 6 packs, 8 packs, and 12 packs.



Thermal Packs In Unit



Standard Pack







Whitehall's Thermalator® Moist Heat Therapy line features stainless steel units which store and heat Thermal Packs at a consistent temperature. Heated Thermal Packs are then applied to affected areas, penetrating heat deep into the muscle to provide therapeutic circulation & pain relief.

Whitehall's Moist Heat Thermal Packs contain bentonite clay filling which retains and evenly distributes heat for up to 30 minutes of therapy. Heavy duty nylon thread keeps them durable while 100% pure cotton keeps them soft & comfortable. Terry covers optional.

Whitehall's Glacier Packs are durable as well as flexible. Constructed with a 13 mil PVC cover, our packs contain flexible filler that never freezes, allowing them to be more effectively molded to the treatment area. Each reusable pack is double seamed for added durability and provides up to 30 minutes of cold therapy treatment.

Whitehall's Glacier Freeze provides convenient and reliable storage for all Glacier Packs. It's designed to keep them properly stored and ready to use. Each unit comes with stainless steel shelves for 12 standard Glacier Packs and a molded compartment for Ice Cups*, which are used to penetrate cold therapy deep into muscles. Oversized packs are easily stored on the bottom shelf.

Chill out with Little Freeze. This compact chilling unit comes standard with specially designed shelves that store four included Glacier Packs, and a molded compartment for Ice Cups*, providing easy access and removal for quick pain relief. The Little Freeze performs the same big chill as the Glacier Freeze in one small, economical unit.

*Ice Cups not included

30LORS & LOGO



SPORTS EQUIPMENT - cleatPRO®


The **cleatPRO**® features an enclosed tank that steams cleats while they hang on an interior rack.

Steaming time ranges between 2-30 minutes based on athlete's preference. cleatPRO $^{\circ}$ can be used to steam and contour the fit of new unworn cleats, as well as for daily steaming of pre-worn cleats. Patent pending.





WHIRLPOOL ACCESSORIES



AS1, AS2 and AS3 Adjustable Suspension Seats



RS Reclining Seat




ARS Arm Rest Support



MXT-15-OTG Thermostatic Mixing Valve



MXWH Washout Hose Attachment



TTS1, TTS2 Tank Top Seat



HR Head Rest, Adjustable



HydroChlor® Whirlpool Antiseptic



When combined with proper cleaning, HydroChlor[®] whirlpool antiseptic packets quickly and easily disinfect whirlpools and ice baths.


Kills:

Reproductive bacteria: 0.2% one(1) hour Tubercle bacillus: 5% two(2) hours Bacillus: More than ten(10) hours



TWO SHELF UTILITY CART

ADDITIONAL PRODUCTS



BOTTLE WARMERS



PARATHERAPY

Manufactured from recycled materials.



